

FIGURE 1





FIGURE 2

Atty Docket: ABIOS.042A



313 To 1

FIGURE 3




FIG. 4A




FIG. 4B

Stockwell et al.
Appl. No.: Unknown Atty Do.




FIG. 4C

Table 1
System components and tunable parameters

Program	Function performed	
BlastParse.pl	Parse BLAST output	
mark_substitution_heteroplasmy.pl	Mark ambiguous base calls	
extract_SE_consensus.pl	Extract a consensus sequence from a phrap polyphred run, replacing low quality bases with '?', and polyphred rank 1 bases with 'N'.	
compute_coverage.pl	Calculates overall, forward strand, and reverse strand coverage for all bases in the consensus sequence	
count_hv1_deletes.pl	Count deletions in HV1 relative to the rCRS	
border_index.pl	Compute start and end positions of HV1 and HV2 regions.	
fix_mitotype_reporting_range.pl	Output a mtDNA profile based on input variants list and reportable range.	
find_bad_traces_from_blast_report.pl	Reject reads that do not align appropriately to the rCRS	
determineReadTypes.pl	Adds template name, template type, and primer type to phred output files.	
seq2delta_vs	Align mtDNA profile to rCRS and report variants according to nomenclature.	
calculate_coverage_mitotype.pl	Mask mtDNA consensus sequence based on required coverage and calculated coverage	
count_hv1_inserts.pl	Count insertions in HV1 relative to the rCRS	
generate_hv_mask_fasta_files.pl	Extract the HV1 portion and HV2 portion of the mtDNA consensus sequence, based on the computed HV1 and HV2 regions.	
flip_fasta.pl	Invert an X masked sequence output by cross_match	

MITOCHONDRIAL DNA AUTOSCORING SYSTEM Stockwell et al.

Appl. No.: Unknown

Atty Docket: ABIOS.042A

Table 2

	Assessment 1	Assessment 2
Total True Positives: N _{automated} = N _{manual}	16	10
Total True Negatives: A _{automated} = A _{manual}	77,358	110,354
Total False Positives: N _{automated} = A _{manual}	95	807
Total False Negatives: A _{automated} = N _{manual}	6 ^a	14
Total Incorrect	0	0
Sensitivity: TP/(TP+FN)	72.73%	41.67%
Specificity: TN/(TN+FP)	99.88%	99.27%
Positive Predictive Value: TP/(TP+FP)	14.41%	1.22%
Negative Predictive Value: TN/(TN+FN)	99.99%	99.99%
,		

^a Predominate base correctly called

Differences take into account the consensus of the two analysts